

# CappMicrobiology

## SINGLE CHANNEL PIPETTES

CappMicrobiology pipettes are ideal for microbiology tests in the food industry.

This pipette is supplied with three fixed volume knobs at 100µl, 900µl and 1000µl and has a modified tip cone suited for use with a sterile extended length straw. This allows further reach and pipetting of viscous liquids such as milk, cream, yoghurt and cosmetics.

It is an ideal instrument for protocols involving dilution of 1:10.

### Highlights:

- Extended Reach
- Rotation Free Volume Selection
- Pre-set Calibrated Volume Controllers
- Fully Autoclavable
- Customizable Volume Selection
- Fully Autoclavable



## Ordering Information

CappMicrobiology pipettes (for use with sterile straws)

Cat. No.	Description	Inaccuracy %	Imprecision %
C1000-1BS	CappMicrobiology pipette, with 100µl, 900µl and 1000µl volume knobs	3.00/1.60	0.80/0.30
4170807	Pre-sterile straws for CappMicrobiology pipette, 40x25 pcs.		