

Lab Burners Dragon 200



Features

- > Automatic ignition due to piezo effect
- > Max. Temperature: 1300°C
- > Adjustable flame
- > Refillable butane tank

Applications

- > Sterilizing
- > Welding
- > Heating
- > Cooking

Specifications

Model	Dragon 200
Max. temperature (°C)	1300
Fuel tanker capacity (g)	45
Working time (min)	60
Dim. H x W x D (mm)	92x115x79
weight (g)	407
Safety switch	--
Order No.	177200-00

WIGGENS GmbH

Gässlesweg 22-24, 75334 Straubenhardt, Germany
Tel.: 0049 7248 4529088

WIGGENS China

Room 426, Hall A, Office Building M8, No.1 Jiuxianqiao East Road,
Chaoyang District, Beijing 100015, China
Tel: +86 400-809-2068
Fax: +86 400-809-2068-112
Email: info@wiggens.com service@wiggens.com
Website: www.wiggens.com