



Roller Apparatus



Celrol Pro

- > Smooth ramping digitally controlled rotation accurate to 0.01 Rpm
- > Maintenance-free precision brushless motor
- > Bright TFT display for easy operation
- > SMS alarm function
- > Can Be directly connected to the printer

Specifications

Model	Celrol Pro
Display	TFT
Capacity	5~55 Positions
Bottle Speed (110mm Bottle)	0.08~2.5rpm
Diameter of roller bottles	108~121mm
Length of roller bottles	Max. 550mm
Accuracy	0.01rpm
Drive	Belt Driven
Motor Type	Brushless DC
Rotation Direction	Clockwise and Counter-Clockwise
Options	Temperature Alarm and Battery Backup
Remote Communication	RS232
Humidity	80% up to 37°C
Operating Temperature	10 to 40°C
Electrical	100-240 VAC, 50/60 Hz

Celrol Elite

- > Robust touch screen interface that is easy to view from a distance
- > Readily recognizable icons for intuitive navigation
- > Digitally controlled maintenance-free motor accurate to 0.01 Rpm
- > RS-422 Modbus RTU, RS-4485 Modbus RTU, Modbus TCP and TCP/IP
- > SMS alarm function
- > Can Be directly connected to the printer

Specifications

Model	Celrol Elite
Display	SIMATIC HMI Touch Screen
Capacity	8~88 Positions 5~55 Positions
Bottle Speed (110mm Bottle)	0.08~2.5rpm
Diameter of roller bottles	108~121mm
Length of roller bottles	Max. 550mm
Accuracy	0.01rpm
Drive	Belt Driven
Motor Type	Brushless DC
Rotation Direction	Clockwise and Counter-Clockwise
Options	Temperature Alarm and Battery Backup
Remote Communication	RS-422 Modbus RTU, RS-4485 Modbus RTU, Modbus TCP and TCP/IP
Humidity	80% up to 37°C
Operating Temperature	10 to 40°C
Electrical	100-240 VAC, 50/60 Hz